

# The **Handi-Grip** Portable Grab Bar

**Handi-Grip** portable grab bars for assistance anywhere in your home or while traveling. The strong suction adheres to most smooth surfaces. Attach **Handi-Grip** portable grab bars in the bath for extra support when exiting the tub, use as a temporary towel bar or to place the toilet paper in a more convenient location. Attach it directly to your tub for help while bathing. Use **Handi-Grip** portable grab bar near your favorite chair for extra assistance when standing, or on a stairwell for an aid climbing steps.

## Features

- Portable, quickly and easily installed and removed
- Telescoping bar for one or two-handed grip
- Large lever for easy release
- Easily attaches to most smooth surfaces



Because the **Hand-Grip** portable grab bars are not screwed into the surface, they are perfect for renters, traveling or visiting. You can easily attach a **Handi-Grip** portable grab bar to marble, glass or tile, surfaces that may be damaged by drilling holes. Take your grab bar when staying in hotels or as a guest in a friend's home. It quickly attaches to a wall and just as quickly comes off to pack in your suitcase. **Handi-Grip** portable grab bars eliminate hiring a contractor for installation. Use a portable grab bar to help determine the best location for a permanent grab bar without drilling trial holes.

Accessory joint adapter allows the grab bar to attach to inside corners or curves. Use the adapter to position the grab bar to the inside of your tub. The aluminum joint swivels a full 360°. The joint adapter fits between the bar and suction housing with spring loaded clips. It is quickly attached and removed from the grab bar when not needed.

## Components

**Handi-grip** portable grab bar consists of two suction cups with release levers and a telescoping bar. The suction housing and bar are high tensile plastic components, suction pads are natural rubber and all metal parts are stainless steel. Available in three telescoping sizes.



**Grab Bar Specialists, Inc.**

a Division of CR Consulting & Distributing, Inc.  
TOLL FREE - 1-800-707-5033

Handi-Grip

# The **Handi-Grip** Portable Grab Bar



## Installations

Due to the vacuum suction generated by the discs, the grab bar can only be attached to smooth and non-porous surfaces, which are pore- and steam-tight, dry and free of oils and grease. On rough or porous surfaces, the vacuum hold diminishes or does not exist at all. All tiles are not necessarily smooth and steam-tight, test your surface first. Using the grab bar under water or in extreme temperatures (such as a sauna) reduces the suction power. The minimum area for each suction cup for adhesion is 5"x5". The use of **Handi-Grip** portable grab bar is not recommended on dry wall.

SPECIFICATIONS	SMALL	MEDIUM	LARGE	SINGLE HAND
overall length	17 1/4" - 22 1/4"	22" - 27"	25 3/4" - 30 3/4"	12 3/4"
length of grab area	7" - 12"	12 5/8" - 17 1/2"	15 3/4" - 20 5/8"	3 3/4"
length center to center of suction cup	13" - 18"	17 1/2" - 22 1/4"	21 1/2" - 26 1/4"	8 1/2"
diameter of suction cup	4 5/8"	4 5/8"	4 5/8"	4 5/8"
outer diameter of bar	1 1/4"	1 1/4"	1 1/4"	1 1/4"
diameter of inner bar that telescopes	1"	1"	1"	NA
height of bar	3"	3"	3"	3"
height of suction cup housing	1"	1"	1"	1"
weight	1 3/4 lb	2 lb	2 1/4 lb	1 1/2 lb
weight capacity *	130 lb	130 lb	130 lb	130 lb
Product number	R1400521	R1400522	R1400523	R1403999

\* **Handi-Grip** portable grab bars are intended to assist the body's natural movement, example: as an aid for standing up. They are not designed to support your whole body weight. The bar is intended to provide support only for short periods of time.

### Grab Bar Specialists, Inc.

a Division of CR Consulting & Distributing, Inc.  
TOLL FREE: 1-800-707-5033

